


Wayflex

Total Wave

Massage Device

ENGLISH




USER'S GUIDE

- ▲ Read all precautions and instructions before using this equipment.
- ▲ After reading, keep this user's guide in a safe place for future reference.

TABLE OF CONTENTS

Important Precautions.....	3
Display and Power Switch Instructions.....	4
Attachment Demonstrations.....	5-6
Battery & Protection System & Charger	7
Speed Data	8
Error Codes and Solutions	9

YOUR MESSAGE DEVICE



Please remember to respect the local regulations: hand in the non-working electrical equipment to an appropriate waste disposal center. The packaging material is recyclable. Dispose of the packaging in an environmentally friendly manner and make it available for the recyclable material collection-service.



Please be sure that you read all warnings on your product before use. Failure to do so may result in injury. Users assumes the risk of injury and all liability resulting the use of this product.

IMPORTANT PRECAUTIONS

To reduce the risk of serious injury, read all important precautions and instructions in this user's guide before use. Users assumes the risk of injury and all liability resulting the use of this product. It is also user's responsibility to take proper care of your product.


WARNING





- Before beginning any exercise program, consult your physician. This is especially important for persons over age 35 or persons with pre-existing health problems.
 - Use the device only as described in this manual.
 - The device is intended for home use only. Do not use the device in a commercial, rental, or institutional setting.
 - The device is for adults only. Children are not suitable to use the device.
 - Keep children under age 12, disabled persons and pets away from the device. Children should be supervised to ensure that they do not play with the appliance.
 - This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
 - Keep the device indoors, away from moisture and dust. Do not put the device in a garage or covered patio, or near water.
 - Inspect and properly tighten all parts regularly. Replace any worn parts immediately.
 - Over exercising may result in serious injury or death. If you feel faint or you experience pain while exercising, stop immediately and cool down.
 - Always check your unit before each use, and ensure that all nuts and bolts are tight, and the unit is in good working condition.
 - Always turn the power off if the device is not being used.
 - Do not disassemble the device by yourself.
- For any technical issues, please contact our service center.

ATTENTION


- It is the responsibility of the owner to ensure that all users of the device are adequately informed of all precautions.
- Warm up at least 5 minutes before workout and cool down at least 5 minutes afterwards.

DISPLAY INSTRUCTIONS




	Press  to turn on/off the device
	Press + to increase the current speed
	Press - to decrease the current speed

POWER SWITCH




ATTACHMENT DEMONSTRATIONS


Belt attachment




Belt attachment with adjustable belt and belt handle




Curved head attachment




Pointed head attachment



Triangle head attachment



Ball head attachment



ATTACHMENT DEMONSTRATIONS



The ball head attachment works on large muscle groups for deep massage for areas like your shoulders, waist, back, and legs.

The triangle head is great for relieving stiff necks and typical shoulder tension and can help relieve headaches, and also it is perfectly shaped to relieve those muscles down your calves or your Achilles area, to help gain mobility and flexibility.



The curved head molds to your body and features stimulation nodules that help lactic acid break down in your stiff muscles, and stimulate blood flow to help increase circulation.

The pointed head is for targeting those really sore areas that need some extra attention. It also works on all those acupuncture points on your body for relief.



The adjustable belt suits whatever length you need for the whole body.

BATTERY


- It takes around 90 min to charge the battery from 0% to 100%.
- A fully charged battery will be able to operate for 5 hours at speed level 01.
- A fully charged battery will be able to operate for 2 hours at speed level 30.
- Please note the speed level will affect the maximum usage time.

PROTECTION SYSTEM

- The device will shut down upon operating for 15 min.
However, you can restart the device immediately without any issues.

CHARGER

- The charger shows RED light when it is charging the battery.
- The charger shows GREEN light when it is full charged.



SPEED DATA

SPEED	VIBRATION PER SECOND
1	14
2	15
3	16
4	17
5	19
6	20
7	21
8	23
9	24
10	26
11	27
12	29
13	30
14	32
15	33
16	34
17	36
18	37
19	39
20	40
21	42
22	43
23	44
24	46
25	47
26	49
27	50
28	52
29	53
30	55

ERROR CODES AND SOLUTIONS

Error code	Solution
E1	The device run out of the battery. Please charge the battery.
E2	Signal failure between PCB board and battery. Please contact the service center.
E3	Motor failure. Please contact the service center.

CAUTION:

Do not disassemble the product by yourself. For any technical issues, please contact our service center.